



Enhanced Sensitivity 4 channel framing camera



Up to 1 Billion frames per second capture speed

50lp/mm and 36lp/mm system resolution options

1360 x 1024 pixel, 12-bit sensor resolution

4 discrete intensified optical channels

The Specialised Imaging SIM-ES offers twice the light transmission compared with standard SIM cameras. Each of the 4 channels provides high resolution images through a standard filter holder making the SIM-ES ideal for multispectral or colour imaging.

Comprehensive triggering adjustment and a wide range of output signals are controlled using the custom software package which also includes measurement and image enhancement functions

FEATURES

- Dedicated standard sized filter holders for each channel
- ☐ Fully adjustable interframe time to 1ns
- ☐ Fully adjustable exposure down to 3ns
- ☐ Gain adjustment up to 10,000X
- ☐ Adjustable output triggers
- □ Nikon lens mount fitting
- ☐ Ethernet communications



Phosphor screen

Enhanced Sensitivity 4 channel framing camera

P46



MODELS SIMX4-ES SIMD8-ES Number of Channels 4 4 Number of Images 4 8 System resolution 50lp/mm 36lp/mm Gain 10,000 10,000

P43

OPTICAL	
Optics	Single input beam splitting optics channels can be easily fitted with standard 31mm dia 1.1mm thick filters
Lenses	Nikon F-Mount
Internal electro- mechanical iris	f2.8 - f22
Shutter	Electro-mechanical
Distortion	Nominally zero
Channel Registration	Within one pixel after software correction
Intensity Variation	Better than 5% across the image

INTENSIFIER / Image Sensor	ICX285AL (Intensified)
Active CCD Pixel	1360 (H) x 1024 (V)
Pixel Size	6.45 µm (H) x 6.45 µm (V)
Dynamic Range	12 bits
Intensifier	Gen II 18mm High resolution MCP Input window Fused Silica Output window Fibre Optic Photocathode S25, others available on request Phosphor Screens: SIMX4-ES: P43 SIMD8-ES: P46 Gen III Intensifiers available on request

MECHANICAL	
Dimensions in cm (LxWxH)	53.3 x 37.6 x 25.5 (without lens, not including periscope +4cm)
Mount	3/8-16 UNC Female
Weight	16Kg approx. (without lens)

TIMING PARAMETERS		
System Clock	1GHz quartz crystal controlled	
Exposure Mode (each image)	Single exposure or multiple exposures (Max. 8) per channel	
Exposure Time	3ns - 10ms in 1ns steps independently variable	
Interframe Time	Ons - 20ms in 1ns steps independently variable	
Delay to 1st exposure	55ns to 10ms in 1ns steps, independently variable	
Flash Outputs	5ns - 1ms in 1ns steps independently variable	
Framing rates	up to 1 Billion fps	

INPUT / OUTPUT SIGNALS		
Trigger 1	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination	
Trigger 2	Electrical signal (BNC connector) Threshold variable from ± 25V Positive or Negative polarity, Make/Break 50Ω or 1KΩ termination	
Timing Monitor Pulses	Pulse width (min. 3ns) and position user programmable TTL into 50Ω	
Flash Trigger Outputs	Pulse width (min. 5ns) and position user programmable TTL into 50Ω	
Camera control	Data and command transfer via 100Mbps ethernet cable length 10m (standard), other lengths up to 100m available 100FX fibre optic ethernet link (up to 2Km) – optional	
Software	Custom software compatible with Microsoft Windows Operating Systems for camera control, image data archiving in various file formats.	
Power Requirements	100-240V AC 2A, 50-60Hz	

ENVIRONMENTAL	
Storage temperature	-10°C to +50°C
Operating temperature	-5°C to +40°C
Humidity	10 - 90% RH non condensing
Vibration shock	10 - 40 Hz Max. 10g in any direction
EMC	Meets all EC harmonized standards

UK (Head Office / Factory)

6 Harvington Park, Pitstone Green Business Park Pitstone. LU7 9GX England

Tel +44 (0) 1442 827728

USA

Specialised Imaging Inc. 40935 County Center Dr. Suite D Temecula, CA 92591, USA

Tel +1 951-296-6406

GERMANY

Hauptstr. 10, 82275 Emmering Germany

Tel +49 8141 666 89 50







FM 87429

ISO9001:2015

specialised-imaging.com

info@specialised-imaging.com